

SSI3 Shaking Incubator



About the SSI3

The SSI3 has a transparent hood that lifts up via a hydraulic system, so it can function in tight places. All of our shaking incubators feature stainless steel interiors which provide excellent durability and an easy-to-clean surface. Each unit has an easy-to-read LED display. The rotation platform is included with each unit and is self-centering for easy installation. The SSI3 includes a convenient, user adjustable counter-balance that provides optimal load flexibility and compensates for off-center loads.

The SSI3 features a small footprint and slim hood design to maximize space and functionality.



Model Name	SSI3 / SSI3-2
Capacity	3.3 cu.ft. (93 liters)
Interior Dimensions (W x D x H)	19.0" x 20.0" x 17.7" 482 x 508 x 450 mm
Exterior Dimensions (W x D x H)	21.0" x 24.0" x 27.0" 533 X 610 X 686 mm
Shaking Speed	30RPM to 400RPM
Temperature Data	
Temperature Range	Ambient +8°C to 60°C
Temperature Uniformity	+/-0.5°C at 30°C



Precise Temperature Control - Superior Uniformity

- Independent Over Temperature Thermostat
- Over Temperature Alarm
- Orbital-Shaking Speed Alarm
- Digital Timer
- Adjustable Orbit

**2 Year
Limited Warranty!**

Accessories

9530526	Flask Holder 500 mL
9530528	Flask Holder 25 mL
9530529	Flask Holder 50 mL
9530530	Flask Holder 125 mL
9530531	Flask Holder 250 mL
9530532	Flask Holder 1000 mL
9530551	Flask Holder 2000 mL
9530553	Flask Holder 2800 mL
9530554	Flask Holder 4000 mL

Accessories

9530555	Flask Holder 6000 mL
9530560	Flask Holder 5000 mL
9751177	10-13 mL Test Tube Shaker (100 tubes)
9751178	14-16 mL Test Tube Shaker (80 tubes)
9751179	18-20 mL Test Tube Shaker (60 tubes)
9751180	22-25 mL Test Tube Shaker (36 tubes)
9751181	50 mL Test Tube Shaker (29 tubes)
9900556	Microtiter Plate Holder (2 per set)
9900687	Flask Holder Bundle (10x25 mL, 10x50 mL, 5x125 mL, 5x250 mL, 2x500 mL, 1x1 L, 1x2 L, 1x4 L)